

# Reusable Digital Forehead Thermometer Products

## Instructions for Use

### How to Read

1. Apply to center of forehead just above the eyebrows and hold in place.
2. Read thermometer when the boxes stop changing color (generally 15 seconds), and while it is still on the forehead.
3. Green indicates the correct temperature. If green does not appear, temperature is between blue and tan.


### Intended Use


Measure human forehead surface temperature and monitor temperature trends

### Accuracy

$\pm 1^{\circ}\text{F}$  ( $0.6^{\circ}\text{C}$ )

### Storage & Handling

 Store at 60 to 104°F (15.6 to 40°C)

 Store at 15 to 95% RH

- Optimal storage at 60 to 80°F (15.6°C to 26.7°C) in 40 to 60% RH

### Precautions

- Environmental temperatures might affect the performance of this product
- Product should be used indoors at normal room temperature



Avoid exposure to UV light



- If cleaning is required, use soap and water or rubbing alcohol



Do not soak

Liquid Crystal Forehead Temperature Trend Monitors are designed to be unbreakable noninvasive, quick reading devices that reliably indicate trends in core body temperature by measuring the surface temperature of the forehead. Liquid Crystal Forehead Temperature Trend Monitors indirectly measure core body temperature by measuring the forehead temperature, whereas clinical thermometers (oral/rectal) directly measure core body temperature. Unusual Liquid Crystal Forehead Temperature Trend Monitor readings should be verified by a clinical thermometer.

For external use only - Conforms to ASTM E-1061-01 - Not made with natural rubber latex

LCR Hallcrest, 1911 Pickwick Lane, Glenview IL 60026 USA  
LCR Hallcrest Ltd. Dock Road, Connah's Quay, UK CH5 4DS

*Hallcrest*